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### **Genomic DNA isolation procedure for FFPE samples**

1. Cut 5-10 sections of up to 10 µm thick and place them in a 2 ml microcentrifuge tube (provided).
2. Add 500 µl of **PRB solution** and vortex vigorously for 10 seconds.
3. Incubate at 56°C for 5 minutes.
4. Centrifuge at 11000g (approx. 10000rpm) for 1 minute at room temperature.
5. Discard the supernatant with care and repeat steps 2-5.
6. Add 500 µl of 100% ethanol and incubate at room temperature for 10 minutes.
7. Centrifuge at 11000g for 5 minutes and discard the supernatant.
8. Repeat steps 6-7 one more time.
9. Dry the sample by incubating the tube opened at 37°C for 5 minutes.
10. Add 300 µl of **Solution A** and vortex gently.  
**Note:** *Mix gently before addition if Solution A was stored at 4°C*
11. Incubate at 56°C for 60 minutes or until the sample has been completely lysed. Use of a horizontal shaker (at 100-150rpm) is optional but preferable. An overnight incubation at 56°C is recommended to ensure complete lysis.
12. Incubate at 90°C for 60 minutes. If using the same heating block, leave the sample at room temperature after the 56°C incubation until the heating block has reached 90°C.  
**Optional RNase treatment:** *heat sample at 70°C for 5 minutes. Add RNase (not provided) 100 µg/ml (final concentration) and incubate at 37°C for 5 minutes. Finally, continue to step 13.*
13. Cool the samples (from step 12) by incubating them approximately 5 minutes on ice.
14. Add 100µl of **Solution B** and vortex gently.
15. Centrifuge at 17000g (approx. 13500rpm) for 10 minutes at room temperature.
16. Carefully pipet the supernatant into a sterile 1.5 ml eppendorf type tube, discarding the pellet.
17. Add 40µl of **Solution C** and 400µl of **Solution D**.
18. Shake slightly by inverting the tube several times until a homogenous solution is observed
19. Incubate the tube at room temperature for 10 minutes in a vertical position.
20. Centrifuge at 11200g (12000rpm) for 5 minutes at room temperature.  
**Note:** *Usually, a little pellet forms.*
21. Discard the supernatant with care.
22. Add 500µl of **Solution E**.
23. Centrifuge at 11200g (12000rpm) for 5 minutes at room temperature.  
**Optional:** *An additional wash with 70% ethanol can be done to avoid salt excess in the final sample.*
24. Discard the supernatant with care and place the tube (containing the pellet) open and upside down, over a filter paper, for 5-10 minutes.
25. Add 20-50µl of **Solution F** and pipet up and down carefully to re-suspend the genomic DNA.  
**Note:** *Nuclease-free water can be used but is not recommended for long time storage.*
26. Optional, incubate the tube at 37°C for 30 minutes to help DNA solubilisation.
27. Use immediately or store at 4°C (if to be used during the next 48 hours) or at -20°C (for longer storage).